

MD Biosciences

Division of Morwell Diagnostics GmbH
Gewerbestr. 9 8132 Egg b. Zürich Switzerland
Tel: +41-1 986 26 28 ♦ Fax: +41-1 986 26 30
For USA/ Canada : (651) 641 1770 ♦ Fax: (651) 641 1773
For Israel Toll Free: 1-800-200 MDB ♦ Fax: (8) 940 9973
For United Kingdom (UK) : 0800 0197 200
e-mail: info@mdbiosciences.com website: www.mdbiosciences.com

Rabbit Polyclonal Antibody to:
Human Procollagenase-2
(MMP-13)

Product Code: 2028019

INTRODUCTION

The polyclonal antibody was raised against recombinant human procollagenase-3 expressed in insect cells. Protein A-purified rabbit anti-human procollagenase-3 polyclonal IgG

PRESENTATION

Purified antibody in 100 µl phosphate buffered saline (20 mM Na-phosphate, pH 7.5, 150 mM NaCl) containing 0.01 % thimerosal.

SPECIFICITY

Based on Western blotting, the antibody reacts with human procollagenase-3. It does not react with human gelatinase A, gelatinase B, MT1-MMP, and the catalytic domains of MT2-MMP, MT3-MMP, MT4-MMP and MT5-MMP.

APPLICATION

♦ Western-Blotting with ECL-Detection
Suggested dilution range: 1:2,500 – 1:8,000

♦ ELISA
Suggested dilution range: 1:x,000 – 1:x,000

STORAGE

The antibody is stable until the expiry date given on the label if stored at -70 °C. Repeated freezing and thawing should be avoided.

EXPIRY DATE

see vial

For Research Use Only