

MD Biosciences

Division of Morwell Diagnostics GmbH
Gewerbestr. 9 8132 Egg b. Zürich Switzerland
Tel: +41-1 986 26 28 Fax: +41-1 986 26 30

For USA/ Canada : (651) 641 1770 Fax: (651) 641 1773

For Israel Toll Free: 1-800-200 MDB Fax: (8) 940 9973

For United Kingdom (UK) : 0800 0197 200

e-mail: info@mdbiosciences.com website: www.mdbiosciences.com

Rat Monoclonal Antibody to:

MOUSE T1/ST2

Purified Antibody

Product Code: 101001

CLONE	DJ8
Ig SUBCLASS	IgG1 (light chain not isotyped)
PRESENTATION	The antibody is supplied in PBS, with 0.1% sodium azide as preservative and purified over protein G-sepharose. The characteristics of each lot are tested by FACS analysis with bone marrow derived mast cells.
CONCENTRATION	See vial for actual concentration
SPECIFICITY	This clone recognizes the membrane anchored murine T1M protein on the surface of T helper 2 cells and mast cells. T1M appears on fetal blood derived mast cell progenitors before they express the Fc γ RI, on IL-3-dependent bone marrow derived mast cells, and on mature peritoneal mast cells. The antibody detects T1S protein consisting only of the extracellular portion of the protein, which is secreted from growth factor and proinflammatory cytokine stimulated murine fibroblasts.
APPLICATION	This antibody allows the identification and purification of murine T helper 2 cells and all forms of murine mast cells.
SUGGESTED USAGE	Flow cytometry, immunoprecipitation
STORAGE	Store at +2 ^o - +8 ^o C
EXPIRY DATE	See label on vial
WARNING	This reagent contains 0.1% sodium azide as preservative. Due to potential hazards arising from the build up of this material in pipes, spent reagent should be disposed of with liberal volumes of water
IMMUNISING ANTIGEN	Eukariotically expressed fusion protein of mouse T1 ectodomain and human immunoglobulin Fc domain
REFERENCES	<ol style="list-style-type: none">1. Moritz et al., 1998. Expression analysis of the soluble and membrane -associated forms of the interleukin-1 receptor-related T1 protein in primary mast cells and fibroblasts <i>Hybridoma</i> 17(2): 107-1162. Moritz et al., 1998. The interleukin-1 receptor-related T1 antigen is expressed on immature and mature mast cells and on fetal blood mast cell progenitors <i>J. Immunol.</i> 161 (9) 4866- 48743. Max Löhning, Arne Stroehmann et al, 1998. T1/ST2 is preferentially expressed on murine Th2 cells, independent of interleukin 4, interleukin 5 and interleukin 10, and important for Th2 effector function <i>Immunology</i> 95 (6): 6930-69354. Damo Xu et al., 1998. Selective Expression of a Stable Cell Surface Molecule on Type 2 but Not Type 1 Helper T Cells <i>J. Exp. Med.</i> 187(5): 787-794

For research use only