

## MD Biosciences

Division of Morwell Diagnostics GmbH  
Gewerbestr. 9 8132 Egg b. Zürich Switzerland  
Tel: +41-1 986 26 28 Fax: +41-1 986 26 30  
For USA/ Canada : (651) 641 1770 Fax: (651) 641 1773  
For Israel Toll Free: 1-800-200 MDB Fax: (8) 940 9973  
For United Kingdom (UK) : 0800 0197 200  
e-mail: [info@mdbiosciences.com](mailto:info@mdbiosciences.com) website: [www.mdbiosciences.com](http://www.mdbiosciences.com)

Mouse Monoclonal Antibody to:  
**Human ST2L**  
Biotinylated

Product Code: 101002B

<b>CLONE</b>	B4E6
<b>Ig SUBCLASS</b>	IgG1
<b>PRESENTATION</b>	The biotinylated antibody is supplied in PBS and purified over protein G-sepharose from concentrated tissue culture supernatant. This antibody is supplied with 0.1% sodium azide as preservative.
<b>CONCENTRATION</b>	See vial for actual concentration
<b>SPECIFICITY</b>	This clone recognizes hST2L specifically on the cell surface of human type 2 (Th2, Tc2, NK2 and NKT2) but not type 1 lymphocytes.
<b>APPLICATION</b>	This antibody allows the identification and purification of human type 2 lymphocyte subsets.
<b>SUGGESTED USAGE</b>	Flow cytometry, FACS, cell sorting, immunohistochemistry with hybridoma culture supernatant.
<b>STORAGE</b>	Store at +2° - +8° C. Protect from light.
<b>EXPIRY DATE</b>	See label on vial
<b>REFERENCES</b>	W.L. Chan et al. 2001. Human IL-18 Receptor and ST2L Are Stable and Selective Markers for the Respective Type 1 and Type 2 Circulating Lymphocytes. J Immunol 2001 167: 1238-1244

*For Research Use Only*