

## MD Biosciences

Division of Morwell Diagnostics GmbH  
Gewerbestr. 9 8132 Egg b. Zürich Switzerland  
Tel: +41-1 986 26 28 Fax: +41-1 986 26 30  
For USA/ Canada : (651) 641 1770 Fax: (651) 641 1773  
For Israel Toll Free: 1-800-200 MDB Fax: (8) 940 9973  
For United Kingdom (UK) : 0800 0197 200  
e-mail: [info@mdbiosciences.com](mailto:info@mdbiosciences.com) website: [www.mdbiosciences.com](http://www.mdbiosciences.com)

Rabbit Polyclonal Antibody to:  
**Human Interleukin- 18 Receptor**  
Neat serum

Product Code: 201006

<b>PRESENTATION</b>	0.5 ml vial, supplied with 0.5 ml of pre-immune matched normal rabbit serum as control. No sodium azide.
<b>SPECIFICITY</b>	Human Interleukin 18 Receptor, a cell surface marker present specifically on Th1/Tc1 cells.
<b>APPLICATION</b>	Use as a specific Th1/Tc1 marker.
<b>SUGGESTED USAGE</b>	Suggested usage at 1/25 and detection with biotinylated goat anti-rabbit IgG followed by fluorochrome conjugated streptavidin for Flow cytometry (FACS). Also suitable for cell sorting and immunohistochemistry using fluorochrome detection.
<b>IMMUNOGEN</b>	A keyhole limpet haemocyanin-conjugated peptide to hIL-18R extracellular domain.
<b>STORAGE</b>	Store at +2 <sup>o</sup> - +8 <sup>o</sup> C for short-term use, at -20 <sup>o</sup> C for extended storage and aliquot to avoid repeated freezing and thawing. Since this preparation contains no preservative, it should be handled under aseptic conditions.
<b>EXPIRY DATE</b>	see vial
<b>REFERENCES</b>	W.L. Chan et al. 2001. Human IL-18 Receptor and ST2L Are Stable and Selective Markers for the Respective Type 1 and Type 2 Circulating Lymphocytes. J Immunol 2001 167: 1238-1244

*For Research Use Only*