

MD Biosciences

Division of Morwell Diagnostics GmbH
Gewerbestr. 9 8132 Egg b. Zürich Switzerland
Tel: +41-1 986 26 28 ♦ Fax: +41-1 986 26 30
For USA/ Canada : (651) 641 1770 ♦ Fax: (651) 641 1773
For Israel Toll Free: 1-800-200 MDB ♦ Fax: (8) 940 9973
For United Kingdom (UK) : 0800 0197 200
e-mail: info@mdbiosciences.com website: www.mdbiosciences.com

RABBIT ANTI-MOUSE COLLAGEN TYPE I, IgG FRACTION

(For research use only)

Catalog N°	MD 20151
Description	Rabbit polyclonal antibodies against Collagen Type I extracted and purified from mouse skin. The IgG fraction is obtained by ion exchange chromatography.
Format	IgG fraction, freeze-dried (0.5 ml). Reconstitute with 0.5 ml di H ₂ O and store aliquots at -20°C.
Stability ;	<u>Unreconstituted</u> : 24 months at +2°/8°C or longer at -20°C avoid freeze-thaw cycles by aliquoting before storage at -20°C and use. <u>Reconstituted</u> : 6 months (mini) at -20°C.
Applications	IFA, RIA, ELISA (tested).
Specificity (% ; at a 1:200 RIA dilution)	Mouse Collagen, Type I : 100 Mouse Collagen, Types II, IV : < 0.1 Mouse Collagen, Type III : < 1.0 Human, Chicken, Rat Collagens, Type I : < 0.1
Use	For immunostaining of extra or intracellular components in light microscopy. For the quantitation of Collagen Type I in biological fluids or in culture supernatants.
Working dilutions	In indirect IFA with fluoresceine anti-rabbit IgG conjugate, use at the following dilution : - 1:40 on frozen mouse skin and liver tissues.